



POLYISOPRENE

SURGICAL GLOVES

Best Performance
Non-latex Synthetic



Medispo® PI Features:



Latex-free: Super soft polyisoprene formulation that is also tough to protect you and patients from pathogens. Exceeds all strength requirements according to USA and EU International standards.



Skin-friendly formulation: Our unique chlorine-free formulation production processes combined with the absence of chemical accelerators known to cause allergic responses mean skin irritations are minimised or eliminated. Tested according to ISO EN 10993.



Easy donning: Skin friendly polymer coated for easy donning even with damp hands saving you time.



Consistent calibrated grip: The outer surface of Medispo surgical gloves are specially treated to provide a consistent grip. The grip is monitored and controlled during production and each lot must pass stringent internal grip controls before being released.



Complies with the following Standards and Regulations:

Quality: ISO 9001, EN ISO 13485

Compliance EU: CE Mark, MDR, EN 455 1-4, EN ISO 374 1-5

Compliance USA: 510K, ASTM D3577

Environment: EN ISO 14001:2015



POLYISOPRENE SURGICAL GLOVES

PRODUCT INFORMATION

Product Name / Description	Medispo® Polyisoprene, Sterile, Hand Specific, 280-300mm
Latex / Latex-free/ material	Latex-free, Polyisoprene
Color	White Cream
Cuff type	Beaded cuff
Accelerators used	ZDEC
Texture / Grip	Micro textured, calibrated grip, medium
Sterile / Non-Sterile / Method	Sterile,E-Beam / Gamma
Powder / Powder-free	Powder-free
Design	Anatomical, forward thumb ball, hand specific
AQL value	0.65
International Product Standards	EN 455 1-4, EN ISO 374 1-5:2016, EN ISO 21420:2019, ASTM D3577, ASTM D6978, ASTM F739, ISO 10993, AS/NZS 4179, ISO 10282, JIS T9107
Quality/Environmental Standards	ISO 9001, ISO 13485:2016, ISO 14001
Shelf Life	3 Years
Compliance to regulations	EU MDR 2017/745, USA 510k, MDSAP
Notified body/ Test Laboratories	DEKRA, SATRA, ARDL, BSI
PPE compliance / Type	PPE Regulation (EU) 2016/425, Type B
Biological Evaluation	EN 10993, Primary Skin Irritation Test, Dermal Sensitization Assay, Cytotoxicity Test
Additional tests	Viral Penetration PPE Compliance (ASTM F1671), Chemotherapy Drugs (ASTM D6978), Chemical Residue



POLYISOPRENE SURGICAL GLOVES

PRODUCT SPECIFICATION

Item No.	MPIG-5.5	MPIG-6.0	MPIG-6.5	MPIG-7.0	MPIG-7.5	MPIG-8.0	MPIG-8.5	MPIG-9.0
Size	5.5#	6.0#	6.5#	7.0#	7.5#	8.0#	8.5#	9.0#
Length(±5mm)	280	280	280	285	285	290	290	300
Thickness(±0.03mm)	Fingertip:0.22		Palm:0.19			Cuff:0.16		

PACKING INFORMATION

Pouch material	Long or square plastic peel pouch		
Packing Quantity	Pouch	Inner box	Shipping Carton
	1pair/pouch	50pairs/box	200pairs/carton

EUROPEAN REPRESENTATION

EC Rep & Registration No.	HBM Medical, Coliemore House, 3 Coliemore Road, Dalkey, Co Dublin, Ireland: 653382
EC Rep Contact Details	INFO@4SHBM.COM, Telephone +353 1 582 2308 , www.4shbm.com

Återförsäljare:

MEDI PLAST

Mediplast AB
Box 1004, 212 10 Malmö
T 020-78 80 35
info.se@mediplast.com
www.mediplast.com

Forhandler:

MEDI PLAST

Mediplast A/S
Marielundvej 46E, 2730 Herlev
T 43 44 40 00
info.dk@mediplast.com
www.mediplast.com

Forhandler:

MEDI PLAST

Mediplast AS
Tollbugata 115, 3041 Drammen
T 32 88 11 00
mediplast@mediplast.no
www.mediplast.com

Jälleenmyyjä:

MEDI PLAST

Mediplast Fenno Oy
PL 153, 01531 Vantaa
Puh: 09 276 360
asiakaspalvelu@mediplast.com
www.mediplast.com