

Ameda® Pearl™ Hospital Grade Breast Pump



Powerful Performance in the Hospital and Home.

Combines powerful performance inside with a protective shell on the outside. Pearl is designed to build milk supply amongst pump dependent, late pre-term and term moms. Pearl's modern design combines the durability needed for the multi-user/hospital environment with lifestyle features moms expect. Pearl will serve the needs of hospitals and is the ideal multi-user product for milk banks and rental stations. Pearl is perfectly positioned to support continuity of care in the Hospital and Home.


Ameda Pearl – Proven Hospital Grade Pumping

Clinical study in process.

- **Proven:** Over 90% of moms surveyed were highly satisfied with Pearl
- **Powerful:** Hospital Grade / Multi-User Performance with two modes of pumping
- **Portable:** Lightweight with rechargeable battery and integrated handle
- **Smart:** Vacusense™ Technology ensures consistent vacuum performance in varying altitudes or changing atmospheric conditions
- **Comfortable:** Fully customizable expression mode with 72 possible combinations of speed and suction settings
- **Durable:** Protective housing designed for easy cleaning between users; 3-year limited warranty
- **Whisper Quiet:** Discreet pumping anywhere
- **Safe:** Compatible with Ameda's Universal HygieniKit™, the only FDA-cleared closed system pump kit proven to protect against contaminants like bacteria, viruses and mold growth

Ameda Pearl – Modern Design

Hospital Grade/Multi-User Proven Durable Performance



Ameda Pearl Technical Specifications

Suction Range	Stimulation: 50 - 150mmHg; 6 user-selectable settings Expression: 30 - 250mmHg; 12 user-selectable settings
Cycles per Minute (CPM)	Stimulation: 80 - 120 CPM; 3 preprogrammed settings Expression: 30 - 48 CPM; 6 user-selectable settings
Power Supply	Medical AC Adapter – Output: 15VDC, 2A, 30W Input: 100-240VAC, 50/60Hz, 1A 2600mAh Lithium-ion rechargeable battery; 3 hours per charge
Unit Weight	1.6kg
Unit Dimensions	20 x 18 x 15cm (L x W x H)
Preprogrammed pumping	Mode auto switches from stimulation to expression at 2 mins Auto Power off at 60 mins Vacusense Technology
Warranty	3-year limited



For more information
about Ameda Pearl



Pearl Pump Options and Accessories

Product	SKU	UPC
Pearl Breast Pump, AC power adapter, North America prong adapter, EN/FR/SP IFU, Quick Start Card	201W01	810725025864
Pearl Breast Pump, AC power adapter, EU & UK prong adapters; International IFU	201U01	810725023624
Trolley	201H01	810725025987
Hard Shell Case	501M03	810725025932
Multi-User Breast Pump Soft Carrying Bag – Detachable shoulder strap and insulated cooling compartment	501M02	810725025611
Universal Diaphragm 2-Pack – 2 Universal Diaphragms	000X02	810725024867
AC Power Adapter North America/AU/Philippines	201X05	810725026045
AC Power Adapter EU/UK	201X06	810725026052
Ameda Pearl Tubing Port Cover	201X08	810725024324
Ameda Pearl Replacement Rubber Feet (set of 4)	201X03	810725026021
Ameda Pearl Replacement Screen Overlay Panel	201X04	810725026038
Ameda Pearl AU Prong Adapter	201X01	810725025895



Ameda Pearl Trolley
201H01




Ameda Pearl Replacement Screen
Overlay Panel
201X04



Ameda Pearl Hard Shell
Transport & Storage Case
501M03



Ameda Pearl AC Power Adapter
EU/UK/Asia
201X06



Ameda Pearl Replacement
Rubber Feet
201X03



Ameda Pearl Tubing
Port Cover
201X08



Ameda Pearl AC Power Adapter
North America/AU/Philippines
201X05



For more info about
Ameda's patented
HygieniKits