


FLERGANGS RESPIRATORSLANGER TIL VOKSNE

- Respiratorslange 19 mm I.D. til voksne
- Flergangs
- Vulkaniserede 22 mm muffer i begge ender
- Materiale: Silikone



Respiratorslange til voksne

Æske = mindste kvantitet/bestilling

Varenummer	Beskrivelse	Længde	Æske/Trsp.karton
6622030	Respiratorslange, flergangs, ID 22mm	30 cm	1/10
6622040	Respiratorslange, flergangs, ID 22mm	40 cm	1/10
6622060	Respiratorslange, flergangs, ID 22mm	60 cm	1/10
6622080	Respiratorslange, flergangs, ID 22mm	80 cm	1/10
6622090	Respiratorslange, flergangs, ID 22mm	90 cm	1/10
6622110	Respiratorslange, flergangs, ID 22mm	110 cm	1/10
6622130	Respiratorslange, flergangs, ID 22mm	130 cm	1/10
6622150	Respiratorslange, flergangs, ID 22mm	150 cm	1/10
6622180	Respiratorslange, flergangs, ID 22mm	180 cm	1/10
6622200	Respiratorslange, flergangs, ID 22mm	200 cm	1/10

VBM

Medioplast A/S
Marielundvej 46E, 2730 Herlev
Tel 43 44 40 00
info.dk@mediplast.com

[mediplast.com](https://www.mediplast.com)