



Forhandler:

**MEDI PLAST**

Medioplast A/S

Marielundvej 46E, 2730 Herlev

T 43 44 40 00

info.dk@mediplast.com

www.mediplast.com



## Technical Sheet

## Surgical skin markers

### Gentian-free ink

#### Reference(s)

**SM01** – Regular tip skin marker with ruler  
**SM03** – Regular fine tip skin marker with ruler  
**SMD02** – Dual fine/regular tip skin marker with ruler

#### Description

Sterile surgical skin markers intended for marking on intact skin. Multiple tips are available depending on the precision required. To ensure accurate measurements, skin marker is provided with a 15 cm ruler, and a printed ruler on the marker body.

	Fine tip	Regular tip
<b>Tip material</b>	Polyacetal resin	Acrylonitrile
<b>Body material</b>	Polypropylene	
<b>Ink</b>	Purple medical-grade ink*	

\*Be aware that ink is free of Gentian violet dye: Gentian violet dye is classified as possibly carcinogenic on basis of sufficient evidence of carcinogenicity in experimental animals to humans by WHO organization.

#### Pack factor

**SM01:** 50 units per carton  
**SM03:** 50 units per carton  
**SMD02:** 100 units per carton

#### Packaging

Individually packed in a FFS packaging

#### Legal Manufacturer

Medline International France SAS - EN ISO 13485 certified

#### Regulation

In compliance with Medical Devices Regulation (EU) 2017/745 and UK MDR  
MD Class: I sterile

- Please refer to complete labeling and instruction for use (if applicable) before using the product.



Latex free



2797



0086



© 2022 Medline Industries LP or one of its subsidiaries. All rights reserved.  
TDS\_SkinMarkersMedline\_EN02