


Sterion® Vessel Loops

- Used in occluding, retracting, and identifying arteries, veins, tendons and nerves in surgical procedures
- Extruded from medical-grade silicone
- Radiopaque
- Intended for single-use, and not made with natural rubber latex
- 2 Loops/Pkg; Ward 10 pcs / Trsp- 10 pcs



Art. No.	Description	Colour	Cross Section	Length
01-1-001	Mini	Red	1.5 x 0.80 mm	406 mm
01-1-002	Mini	Blue	1.5 x 0.80 mm	406 mm
01-1-003	Mini	Yellow	1.5 x 0.80 mm	406 mm
01-1-004	Mini	White	1.5 x 0.80 mm	406 mm
01-1-011	Maxi	Red	2.5 x 1.0 mm	406 mm
01-1-012	Maxi	Blue	2.5 x 1.0 mm	406 mm
01-1-013	Maxi	Yellow	2.5 x 1.0 mm	406 mm
01-1-014	Maxi	White	2.5 x 1.0 mm	406 mm
01-1-021	Supermaxi	White	5 x 1.3 mm	559 mm
01-1-022	Supermaxi	Blue	5 x 1.3 mm	559 mm
01-1-032	Micro	Blue	0.53 x 0.53 mm	406 mm
01-1-301	Mini	Red	1.5 x 0.80 mm	762 mm
01-1-302	Mini	Blue	1.5 x 0.80 mm	762 mm
01-1-303	Mini	Yellow	1.5 x 0.80 mm	762 mm
01-1-304	Mini	White	1.5 x 0.80 mm	762 mm
01-1-311	Maxi	Red	2.5 x 1.0 mm	762 mm
01-1-312	Maxi	Blue	2.5 x 1.0 mm	762 mm
01-1-313	Maxi	Yellow	2.5 x 1.0 mm	762 mm
01-1-314	Maxi	White	2.5 x 1.0 mm	762 mm

MDD class IIa